## Celle / Cells

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*The migration test by cell is the most suitable method for inhomogeneous packaging such as multilayer materials (sheets or films) and allows to test with the food simulant only the surface which is intended to come into contact with food (single-side).*

*Migration tests by cell are also suitable for large containers for which the test by filling is not practicable because of their large capacity.*

*In these cases, representative specimens of the sample can be tested.*

*Migration tests by cell have the advantage of readily reproducible sample geometry.*

*Preparing the cells supplied by Sepack Lab by placing the specimens to be tested on both sides of the AISI 304 / X5CrNi18-10 stainless steel ring, a surface to volume ratio of 2 dm<sup>2</sup> of surface area in contact for about 190 ml of food simulant is obtained.*

*The cells supplied by Sepack Lab are equipped with PTFE discs which, in addition to guaranteeing the seal of the sample on the ring, avoid the possible contact of the simulant with the other components of the cell in case of permeation of the food simulants during the migration test which may cause possible interferences in the test results.*

**Cleaning and correct use of the cells:**

*Before using the cell for the first time, carefully clean all its parts with a degreasing solvent (e.g. acetone) and other detergents in order to remove any residual lubricants used for its manufacture.*

*The stainless steel ring can be further cleaned with diluted nitric acid.*

*Never use migration cells above the boiling point of the used simulant. For example, the upper temperature limits defined by the standardized methods for overall migrations are as follows:*

- UNI EN 1186-4:2003 "Test methods for overall migration into olive oil by cell": up to 100 °C;
- UNI EN 1186-5:2003 "Test methods for overall migration into aqueous food simulants by cell": up to 70 °C (except for particular cases in which temperatures higher than 70 °C can be reached).