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*The migration tests by total immersion is the most suitable method for mono-material packaging in the form of films and sheets.*

*Migration tests by total immersion are also suitable for large containers for which the test by filling is not practicable because of their large capacity.*

*In these cases, representative specimens of the sample can be tested.*

*The specimens supports supplied by Sepack Lab are able to hold up and keep the specimens separated while ensuring complete contact with the simulant.*

*The specimen supports are made of an AISI 316 / X2CrNiMo17-12-2 stainless steel wire with a diameter of 2 mm and the crosspieces are fixed by AISI 316 stainless steel welding.*

*The AISI 316 stainless steel welding of the joints limits the corrosion due to the aqueous food simulant acetic acid, 3% (w/v), allowing to obtain a dry residue lower than 5 mg/l, requirement for the blank reagent of aqueous food simulants according to the standardized method UNI EN 1186-1:2003 "Guide to the selection of conditions and test methods for overall migration".*

*In order to ensure that the specimens are well separated and that the surfaces are freely exposed to the food simulant during the test period, especially for thin films, it is advisable to insert a piece of fine mesh stainless steel net between the specimens.*

*Fine mesh stainless steel nets are not supplied by Sepack Lab.*

*Cleaning of the specimens supports:*

*Before using it for the first time, carefully clean the stainless steel specimens support in order to remove any residual lubricants used for its manufacture. The use of a degreasing solvent (e.g. acetone), followed by diluted nitric acid has proven to be suitable.*